

Colonoscopic Polypectomy - Summary Clinical Guideline

Reference no.: CG/2023/037

The ERBE unit is configured with 3 polypectomy settings as defined below. Switch on the machine and choose the polypectomy program required.

Program 1 Blue Pedal, Forced Coag 20W, Effect 2

Program 2 Yellow Pedal Endo Q Effect 1, Forced Coag 15W

Program 3 Yellow Pedal Endo Q Effect 2, 20W

Size	Morphology	Technique	ERBE setting
≤2mm	Any	Cold Biopsy Forceps	-
≤6mm	Any	Cold snare	-
Any size >6mm	Pedunculated (not semi-pedunculated)	Hot snare	Program 1 POLYPECTOMY Blue Pedal Forced Coag 20W, Effect 2
7-15mm	Sessile or semi-pedunculated	Saline EMR +/- methylene blue	Program 3 - EMR Yellow pedal (keep pressed) Endo Q Eff 2, 20W
≥15mm	Sessile or semi-pedunculated	EMR using methylene blue/Indigocarmine	Program 3 - EMR Yellow pedal (keep pressed) Endo Q Eff 2, 20W

Right colon polyps - **Program 2** - Endo Q Effect 1, Forced Coag 15W, (or Program 3)